



LARGE MEDICAL GRADE POLYURETHANE TUBING

This extremely kink-free 85 Shore "A" durometer ether based tubing is very flexible and resistant to abrasion and scuffing. This ether based tubing meets USP Class VI testing and is sterilizable by gamma or chemical. Do not autoclave

	Cat. #	Dimensions I.D." x O.D."	Dimensions I.D.mm x O.D.mm
	BB520-63	0.0625 x 0.125	1.58 x 3.18
	BB520-94	0.0937 x 0.189	2.40 x 4.88
	BB520-125	0.125 x 0.250	3.18 x 6.35
	BB520-188	0.188 x 0.313	4.78 x 7.95
	BB520-250	0.250 x 0.375	6.35 x 9.53

GENERAL PURPOSE FDA APPROVED POLYURETHANE

This ether-based polyurethane tubing has a hardness of 82 Shore "A" durometer. This tubing is very flexible and resistant to abrasion and tearing. It is extremely resistant to water, air, chemical and petroleum based products. This ether based tubing is FDA compliant for food contact. It will withstand temperatures of -100°F to 175°F.


Cat. #	Dimensions I.D." x O.D."	Dimensions I.D.mm x O.D.mm	
BB521-63	0.063 x 0.125	1.60 x 3.17	
BB521-125	0.125 x 0.188	3.17 x 4.78	
BB521-188	0.188 x 0.250	4.78 x 6.35	
BB521-250	0.250 x 0.375	6.35 x 9.53	
BB521-313	0.313 x 0.438	7.95 x 11.13	
BB521-375	0.375 x 0.500	9.53 x 12.70	
BB521-500	0.500 x 0.625	12.70 x 15.88	


GENERAL PURPOSE FDA APPROVED POLYURETHANE


This ether-based polyurethane tubing has a hardness of 82 Shore "A" durometer. This tubing is very flexible and resistant to abrasion and tearing. It is extremely resistant to water, air, chemical and petroleum based products. This ether based tubing is FDA compliant for food contact. It will withstand temperatures of -100°F to 175°F.

	Cat. #	Dimensions I.D." x O.D."	Dimensions I.D.mm x O.D.mm
	BB523-050	0.063 x 0.125	1.60 x 3.17
	BB523-063	0.125 x 0.188	3.17 x 4.78
	BB523-125	0.125 x 0.250	3.17 x 6.35
	BB523-188	0.188 x 0.313	4.78 x 7.95
	BB523-250	0.250 x 0.375	6.35 x 9.53
	BB523-313	0.313 x 0.500	7.95 x 12.70
	BB523-375	0.375 x 0.500	9.53 x 12.70
BB523-500	0.500 x 0.750	12.70 x 19.05	


CATHETER CONNECTORS


POLYACETAL THREADED FEMALE LUER LOCK		
SCI Cat. #	Description	Pack of 20
BB9257	1.40mm Through Hole	
BB9258	3.10mm Through Hole	

POLYAMIDE THREADED FEMALE LUER CONNECTOR			
SCI Cat. #	Description	Color	Pack of 20
BB9290	1.25mm Through Hole	White	
BB9291	1.75mm Through Hole	White	
BB9292	2.05mm Through Hole	White	
BB9293	2.40mm Through Hole	Blue	

POLYCARBONATE CATHETER HUB		
SCI Cat. #	Connects Tube to Female Luer	Packs of 20
BB9270	.056 - .058" (4 French)	
BB9271	.069 - .071" (5 French)	
BB9272	.081 - .083" (6 French)	
BB9273	.107 - .109" (8 French)	

POLYAMIDE FEMALE LUER BALLON CATHETER		
SCI Cat. #	Description	Packs of 20



BB9262	1.30mm Through Hole. Threaded end and lock wings.	
BB9263	1.55mm Through Hole. Threaded end and lock wings.	

HDPE STEPPED WINGED CATHETER ADAPTER		
SCI Cat. #	Connects	Packs of 20
BB9260	.110 to .281" O.D. Female Luer	
BB9261	.203 to .484" O.D. Female Luer	

POLYPROPYLENE TUBING CONNECTORS (SMALL)

These small diameter polypropylene tubing connectors have excellent chemical resistance. They can be easily steGamma radiation or autoclaved at 250°F. They are designed to withstand 20 p.s.i.






MALE LUER WITH LOCK RING			
Cat. #	Connects	Price	Packs of 12
BB822B	3/32" I.D.	RMB120.00	
BB822C	1/8" I.D.	RMB130.00	
BB822D	5/32" I.D.	RMB140.00	
BB822E	3/16" I.D.	RMB160.00	
BB822F	1/4" I.D.	RMB160.00	

"T" TUBE TO TUBE CONNECTORS			
SCI Cat. #	Connects	Price	Packs of 12
BB825A	1/16" I.D. 3-Way	RMB160.00	 
BB825B	3/32" I.D. 3-Way	RMB160.00	
BB825C	1/8" I.D. 3-Way	RMB160.00	
BB825D	5/32" I.D. 3-Way	RMB170.00	
BB825E	3/16" I.D. 3-Way	RMB190.00	
BB825F	1/4" I.D. 3-	RMB190.00	

	Way			
--	-----	--	--	--

EXTENDED MALE AND FEMALE COUPLER Pack of 12

Cat. #	Connects	Price	Cat. #	Description
BB818A	1/16" I.D. Tube	RMB180.00	BB816	1 Inch Long Extended Length Coupler
BB818B	3/32" I.D. Tube	RMB180.00		
BB818C	1/8" I.D. Tube	RMB180.00		
BB818D	5/32" I.D. Tube	RMB190.00		
BB818E	3/16" I.D. Tube	RMB190.00		
BB818F	1/4" I.D. Tube	RMB190.00		

DUAL MALE LUER CONNECTOR				FEMALE LUER CONNECTORS		
Cat. #	Connects	Pack of 6 = RMB250.00		Cat. #	Connects	Packs of 12 = RMB150.00
BB819	Two male luer connectors			BB9217	Female coupler	
MALE AND FEMALE LUER						
Cat. #	Connects	Pack of 6 = RMB250.00		Cat. #	Connects	Pack of 12 = RMB150.00
BB820	Connects male to female			BB9218	"T" style luer coupler	
ROTATING MALE COUPLER						
Cat. #	Connects	Pack of 6 = RMB450.00		Cat. #	Connects	Pack of 12 = RMB150.00
BB821	Male luer on both ends to female luer connector			BB9219	"L" style luer coupler	